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Cíle

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```
set stepper using units rev/s
for //control.Node :
  configure port A ✓
  rotate stepper through
    angle (rev) 0.8
    to max ±speed (rev/s) -0.33
    with acceleration (rev/s²) 0.1
  Wait for completion ✓
sleep for 3000 ms
set stepper using units rev/s
for //control.Node :
  configure port A ✓
  rotate stepper through
    angle (rev) 0.8
    to max ±speed (rev/s) 0.33
    with acceleration (rev/s²) 0.1
  Wait for completion ✓
```

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```
set //control.Node Speaker (0-20000) frequency to 100 Hz
set stepper using units rev/s
for //control.Node :
  configure port A ✓
  rotate stepper through
    angle (rev) 0.8
    to max ±speed (rev/s) -0.33
    with acceleration (rev/s²) 0.1
  Wait for completion ✓
turn //control.Node Speaker on: false
sleep for 3000 ms
turn //control.Node Speaker on: true
set stepper using units rev/s
for //control.Node :
  configure port A ✓
  rotate stepper through
    angle (rev) 0.8
    to max ±speed (rev/s) 0.33
    with acceleration (rev/s²) 0.1
  Wait for completion ✓
turn //control.Node Speaker on: false
```

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```

set stepper using units rev/s
for //control.Node :
configure port A ✓
rotate stepper through
  angle (rev) 0.8
  to max ±speed (rev/s) -0.33
  with acceleration (rev/s²) 0.1
Wait for completion
beepWhileBridgelsMoving
sleep for 3000 ms
set stepper using units rev/s
for //control.Node :
configure port A ✓
rotate stepper through
  angle (rev) 0.8
  to max ±speed (rev/s) 0.33
  with acceleration (rev/s²) 0.1
Wait for completion
beepWhileBridgelsMoving

to beepWhileBridgelsMoving
set bridgelsMoving to true
repeat while bridgelsMoving
do
set //control.Node Speaker (0-20000) frequency to 100 Hz
sleep for 500 ms
turn //control.Node Speaker on: false
sleep for 500 ms
set bridgelsMoving to value of Angular Velocity rev/s ≠ 0
    
```

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Telefon (USA) 1-800-772-8700 (Volitelné4)

Telefon (mezinárodní) +1 916 462 8384

Online **[pasco.com/support](https://www.pasco.com/support)**